



## SmartHub

### Quick Start Guide

---

#### Tips for Better Results

- Place your SmartHub close to a window or someplace in your house with good cell phone reception.
- Place your SmartHub within 25 feet of where you will take your measurements with your devices.

#### Setup

1. Plug in your SmartHub.
2. Wait for the SmartHub to make two announcements. First, it will tell you that it is powering up. Second, it will tell you that it is ready to send data to your provider. This can take up to 3 minutes.
3. That's it. Your SmartHub is ready for use. There is no software or apps to install!

Your SmartHub **DOES NOT** use your Wi-Fi.

Your SmartHub **DOES NOT** listen for verbal commands. It can talk, but it does not listen.

If your SmartHub does not announce that it is ready within three minutes, look at the bottom of your SmartHub for the Item number and call 855-477-7000. A friendly technician will assist you.

### What It Does

Your RCP SmartHub collects data from the medical devices assigned to you by your provider. It then forwards that data to the practice. The data is encrypted and fully secure. Your SmartHub does not use your Wi-Fi. It works just like a cell phone to send the data.

### Using Your Devices with Your SmartHub

You will be using one or more of these devices:

- Scale
- Blood Pressure Monitor
- Thermometer
- Pulse Oximeter
- Spirometer
- Glucometer

After each measurement, the SmartHub will announce that it has successfully delivered the data to your provider.

### If You Have Any Issues or Questions

Call Customer Support at 855-477-7000.